

## Product Datasheet

### Rat RBP4 ELISA Kit orb411153

<b>Description</b>	Rat Retinol binding protein 4,RBP-4 ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	39 ng/ml-2500 ng/ml.
<b>Sensitivity</b>	39 ng/ml.
<b>Tested Applications</b>	ELISA
<b>Target</b>	RBP4
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>P04916</b>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm