

Product Datasheet

ErbB-2 antibody orb254372

Description	Mouse monoclonal antibody to ErbB-2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide
Target	ErbB-2
Form/Appearance	PBS (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol
Concentration	1.0 mg/ml
Storage	Aliquot and store at -20°C.
Note	For research use only.
Isotype	IgG
Clonality	Monoclonal
Purity	Affinity purification
Expiration Date	12 months from date of receipt.