

Product Datasheet

CD23 antibody orb254375

Description	Mouse monoclonal antibody to CD23
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide
Target	CD23
Form/Appearance	PBS (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol
Concentration	1.0 mg/ml
Storage	Aliquot and store at -20°C.
Note	For research use only.
Isotype	IgG
Clonality	Monoclonal
Purity	Affinity purification
Expiration Date	12 months from date of receipt.