

Product Datasheet

PIWIL2 antibody orb254851

Description	Rabbit polyclonal antibody to PIWIL2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide
Target	PIWIL2
Form/Appearance	PBS (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol
Concentration	1.0 mg/ml
Storage	Aliquot and store at -20°C.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Expiration Date	12 months from date of receipt.