

Product Datasheet

PPP2CB antibody orb254879

| | |
|----------------------------|---------------------------------------------------------------------------------------------------------------|
| Description | Rabbit polyclonal antibody to PPP2CB |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IHC-P, WB |
| Immunogen | KLH conjugated synthetic peptide |
| Target | PPP2CB |
| Form/Appearance | PBS (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol |
| Concentration | 1.0 mg/ml |
| Storage | Aliquot and store at -20°C. |
| Note | For research use only. |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | Affinity purification |
| Expiration Date | 12 months from date of receipt. |