

## Product Datasheet

### CHAF1B antibody orb255069

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Rabbit polyclonal antibody to CHAF1B  |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Human   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | ICC, IF, WB   |
| <b>Immunogen</b>           | KLH conjugated synthetic peptide  |
| <b>Target</b>              | CHAF1B  |
| <b>Form/Appearance</b>     | PBS (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol |
| <b>Concentration</b>       | 1.0 mg/ml   |
| <b>Storage</b>             | Aliquot and store at -20°C.   |
| <b>Note</b>                | For research use only.  |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Purity</b>              | Affinity purification   |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |