

## Product Datasheet

### Human FGFR1 ELISA Kit orb437001

<b>Description</b>	Human Fibroblast Growth Factor Receptor 1 (FGFR1) Sandwich ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156-10 ng/mL
<b>Sensitivity</b>	0.057 ng/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	FGFR1
<b>Storage</b>	All reagents to be kept according to labels on the vials. TMB Substrate, Wash Buffer (30X concentrate), Stop Solution at 4°C; all others at -20°C.
<b>Note</b>	For research use only.
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	tissue homogenates, cell lysates and other biological fluids.