

Product Datasheet

C21orf33 antibody orb255191

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to C21orf33 |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ICC, IF, WB |
| Immunogen | KLH conjugated synthetic peptide |
| Target | C21orf33 |
| Form/Appearance | PBS (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol |
| Concentration | 1.0 mg/ml |
| Storage | Aliquot and store at -20°C. |
| Note | For research use only. |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | Affinity purification |
| Expiration Date | 12 months from date of receipt. |