

Product Datasheet

TSPAN13 antibody orb255816

Description	Rabbit polyclonal antibody to TSPAN13
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide
Target	TSPAN13
Form/Appearance	PBS (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol
Concentration	2.1 mg/ml
Storage	Aliquot and store at -20°C.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Expiration Date	12 months from date of receipt.