

Product Datasheet

Cannabinoid R 1 antibody orb255852

Description	Rabbit polyclonal antibody to Cannabinoid R 1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide
Target	Cannabinoid R 1
Form/Appearance	PBS (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol
Concentration	0.33 mg/ml
Storage	Aliquot and store at -20°C.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Expiration Date	12 months from date of receipt.