

Product Datasheet

GST-Tag antibody orb256055

Description	Rabbit polyclonal antibody to GST-Tag
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Target	GST-Tag
Form/Appearance	PBS (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol
Storage	Aliquot and store at -20°C.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Expiration Date	12 months from date of receipt.