

Product Datasheet

RFP antibody orb256064

Description	Mouse monoclonal antibody to RFP
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Target	RFP
Form/Appearance	PBS (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol
Storage	Aliquot and store at -20°C.
Note	For research use only.
Isotype	IgG
Clonality	Monoclonal
Purity	Affinity purification
Expiration Date	12 months from date of receipt.