

Product Datasheet

c-myc tag antibody orb256359

Description	Mouse monoclonal antibody to c-myc tag
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, IP, WB
Immunogen	A synthetic peptide corresponding to aa 408-439 from C-terminus of human c-myc.
Target	c-myc tag. Epitope sequence, EQKLISEEDL, located in aa 410-419 of human c-myc protein.
Form/Appearance	PBS with preservative (0.02% Proclin 300)
Storage	Store at 4°C for up to 3 months. For longer term storage aliquot into small volumes and store at -20°C.
Note	For research use only.
Isotype	IgG2C kappa
Clonality	Monoclonal
Purity	Protein A affinity purified
Clone ID	90000000000
Uniprot ID	P01106
Expiration Date	12 months from date of receipt.