

Product Datasheet

His tag antibody orb256370

Description	Human monoclonal antibody to His tag
Species/Host	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Purified his tagged Fab fragment.
Target	His tag. Carboxy-terminal oligo-histidine tags (6x His tags) of fusion proteins.
Form/Appearance	PBS with preservative (0.02% Proclin 300)
Storage	Store at 4°C for up to 3 months. For longer term storage aliquot into small volumes and store at -20°C.
Note	For research use only.
Isotype	IgG4 kappa
Clonality	Monoclonal
Purity	Protein A affinity purified
Clone ID	3D5
Expiration Date	12 months from date of receipt.