

## Product Datasheet

### Fluorescein antibody orb256386

<b>Description</b>	Human monoclonal antibody to Fluorescein
<b>Species/Host</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FACS, WB
<b>Immunogen</b>	Fluorescein I-bovine serum albumin (FII-BSA).
<b>Target</b>	Fluorescein. Fluorescein labelled proteins.
<b>Form/Appearance</b>	PBS with preservative (0.02% Proclin 300)
<b>Storage</b>	Store at 4°C for up to 3 months. For longer term storage aliquot into small volumes and store at -20°C.
<b>Note</b>	For research use only.
<b>Isotype</b>	IgG3 kappa
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Protein A affinity purified
<b>Clone ID</b>	4-4-20(enhanced)
<b>Expiration Date</b>	12 months from date of receipt.