

Product Datasheet

Fluorescein antibody orb256388

Description	Mouse monoclonal antibody to Fluorescein
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	FACS, WB
Immunogen	Fluorescein I-bovine serum albumin (FII-BSA).
Target	Fluorescein. Fluorescein labelled proteins.
Form/Appearance	PBS with preservative (0.02% Proclin 300)
Storage	Store at 4°C for up to 3 months. For longer term storage aliquot into small volumes and store at -20°C.
Note	For research use only.
Isotype	IgG2A kappa
Clonality	Monoclonal
Purity	Protein A affinity purified
Clone ID	4-4-20(enhanced)
Expiration Date	12 months from date of receipt.