

Product Datasheet

hapten 4-hydroxy-3-NP acetyl antibody orb256414

Description	Rat monoclonal antibody to hapten 4-hydroxy-3-NP acetyl
Species/Host	Rat
Conjugation	Unconjugated
Tested Applications	FACS, WB
Immunogen	hapten-protein conjugated NP12-chicken gammaglobulin (CG).
Target	hapten 4-hydroxy-3-NP acetyl. Hapten (4-hydroxy-3-nitrophenyl) acetyl (NP).
Form/Appearance	PBS with preservative (0.02% Proclin 300)
Storage	Store at 4°C for up to 3 months. For longer term storage aliquot into small volumes and store at -20°C.
Note	For research use only.
Isotype	IgG2C kappa
Clonality	Monoclonal
Purity	Protein A affinity purified
Clone ID	B1-8
Expiration Date	12 months from date of receipt.