


Product Datasheet

V5 tag antibody orb256432

Description	Mouse monoclonal antibody to V5 tag
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IP, RIA, WB
Immunogen	Pk (V5) Tag (GKPIPPLLGLDST).
Target	V5 tag. Recognises a 14 amino acid sequence (GKPIPPLLGLDST) derived from the P and V proteins of the paramyxovirus simian virus 5.
Form/Appearance	PBS with preservative (0.02% Proclin 300)
Storage	Store at 4°C for up to 3 months. For longer term storage aliquot into small volumes and store at -20°C.
Note	For research use only.
Isotype	IgG1 kappa
Clonality	Monoclonal
Purity	Protein A affinity purified
Clone ID	SV5-P-K
Expiration Date	12 months from date of receipt.



ELISA of anti-epitope tag antibodies on multi-tag protein