

## Product Datasheet

### V5 tag antibody orb256445

<b>Description</b>	Rat monoclonal antibody to V5 tag
<b>Species/Host</b>	Rat
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, ICC, IF, IHC-Fr, IP, RIA, WB
<b>Immunogen</b>	Pk (V5) Tag (GKPIPPELLGLDST).
<b>Target</b>	V5 tag. Recognises a 14 amino acid sequence (GKPIPPELLGLDST) derived from the P and V proteins of the paramyxovirus simian virus 5.
<b>Form/Appearance</b>	PBS with preservative (0.02% Proclin 300)
<b>Storage</b>	Store at 4°C for up to 3 months. For longer term storage aliquot into small volumes and store at -20°C.
<b>Note</b>	For research use only.
<b>Isotype</b>	IgG1 kappa
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Protein A affinity purified
<b>Clone ID</b>	SV5-P-K
<b>Expiration Date</b>	12 months from date of receipt.