

## Product Datasheet

### Mouse Fab fragment kappa isotype control orb256450

<b>Description</b>	Mouse Fab fragment kappa isotype control
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FACS, WB
<b>Target</b>	Mouse Fab fragment kappa. An antibody of Mouse Fab fragment kappa specificity that maybe used as a negative control.
<b>Form/Appearance</b>	PBS with preservative (0.02% Proclin 300)
<b>Storage</b>	Store at 4°C for up to 3 months. For longer term storage aliquot into small volumes and store at -20°C.
<b>Note</b>	For research use only.
<b>Isotype</b>	Fab fragment kappa
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified by Immobilized Metal Affinity Chromatography
<b>Clone ID</b>	MOPC-21
<b>Expiration Date</b>	12 months from date of receipt.