

Product Datasheet

Mouse IgG3 kappa isotype control orb256461

Description	Mouse IgG3 kappa isotype control
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	FACS, WB
Target	Mouse IgG3 kappa. An antibody of Mouse IgG3 kappa specificity that maybe used as a negative control.
Form/Appearance	PBS with preservative (0.02% Proclin 300)
Storage	Store at 4°C for up to 3 months. For longer term storage aliquot into small volumes and store at -20°C.
Note	For research use only.
Isotype	IgG3 kappa
Clonality	Monoclonal
Purity	Protein A affinity purified
Clone ID	MOPC-21
Expiration Date	12 months from date of receipt.