

Product Datasheet

Beta galactosidase antibody orb256417

| | |
|----------------------------|---|
| Description | Human monoclonal antibody to Beta galactosidase |
| Species/Host | Human |
| Conjugation | Unconjugated |
| Tested Applications | FACS, WB |
| Immunogen | beta-Galactosidase. |
| Target | Beta galactosidase. 13R4 is a hyper-stable single-chain antibody directed against beta-Galactosidase. |
| Form/Appearance | PBS with preservative (0.02% Proclin 300) |
| Storage | Store at 4°C for up to 3 months. For longer term storage aliquot into small volumes and store at -20°C. |
| Note | For research use only. |
| Isotype | IgG1 lambda |
| Clonality | Monoclonal |
| Purity | Protein A affinity purified |
| Clone ID | 13R4 |
| Uniprot ID | P16278 |
| Expiration Date | 12 months from date of receipt. |