

## Product Datasheet

### Beta galactosidase antibody orb256420

<b>Description</b>	Human monoclonal antibody to Beta galactosidase
<b>Species/Host</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FACS, WB
<b>Immunogen</b>	beta-Galactosidase.
<b>Target</b>	Beta galactosidase. 13R4 is a hyper-stable single-chain antibody directed against beta-Galactosidase.
<b>Form/Appearance</b>	PBS with preservative (0.02% Proclin 300)
<b>Storage</b>	Store at 4°C for up to 3 months. For longer term storage aliquot into small volumes and store at -20°C.
<b>Note</b>	For research use only.
<b>Isotype</b>	IgG3 lambda
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Protein A affinity purified
<b>Clone ID</b>	13R4
<b>Uniprot ID</b>	<b>P16278</b>
<b>Expiration Date</b>	12 months from date of receipt.