

Product Datasheet

Beta galactosidase antibody orb256425

Description	Mouse monoclonal antibody to Beta galactosidase
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	FACS, WB
Immunogen	beta-Galactosidase.
Target	Beta galactosidase. 13R4 is a hyper-stable single-chain antibody directed against beta-Galactosidase.
Form/Appearance	PBS with preservative (0.02% Proclin 300)
Storage	Store at 4°C for up to 3 months. For longer term storage aliquot into small volumes and store at -20°C.
Note	For research use only.
Isotype	IgG2B lambda
Clonality	Monoclonal
Purity	Protein A affinity purified
Clone ID	13R4
Uniprot ID	P16278
Expiration Date	12 months from date of receipt.