

## Product Datasheet

### AIM1L antibody orb258047

<b>Description</b>	Rabbit polyclonal antibody to AIM1L.
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC-P, WB
<b>Immunogen</b>	Synthetic peptide
<b>Target</b>	AIM1L. This antibody has specificity with AIM1L
<b>Form/Appearance</b>	PBS with 0.02% Sodium Azide, 2% BSA, PH 7.2
<b>Concentration</b>	0.5 mg/ml
<b>Storage</b>	Store at -20°C, Avoid freeze/thaw cycles
<b>Note</b>	For research use only.
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Antigen Affinity Purified
<b>Uniprot ID</b>	<b>Q8N1P7</b>
<b>Expiration Date</b>	12 months from date of receipt.