

Product Datasheet

ywhag antibody orb258316

| | |
|----------------------------|---|
| Description | Rabbit polyclonal antibody to ywhag. |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IHC-P, WB |
| Immunogen | Synthetic peptide |
| Target | ywhag. This antibody has specificity with ywhag |
| Form/Appearance | PBS with 0.02% Sodium Azide, 2% BSA, PH 7.2 |
| Concentration | 0.5 mg/ml |
| Storage | Store at -20°C, Avoid freeze/thaw cycles |
| Note | For research use only. |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | Antigen Affinity Purified |
| Expiration Date | 12 months from date of receipt. |