

Product Datasheet

ZKSCAN5 antibody orb258333

Description	Rabbit polyclonal antibody to ZKSCAN5.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic peptide
Target	ZKSCAN5. This antibody has specificity with ZKSCAN5
Form/Appearance	PBS with 0.02% Sodium Azide, 2% BSA, PH 7.2
Concentration	0.5 mg/ml
Storage	Store at -20°C, Avoid freeze/thaw cycles
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Antigen Affinity Purified
Expiration Date	12 months from date of receipt.