

## Product Datasheet

### RAI14 antibody orb258361

<b>Description</b>	Rabbit polyclonal antibody to RAI14.
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	Synthetic peptide
<b>Target</b>	RAI14. This antibody has specificity with RAI14
<b>Form/Appearance</b>	PBS with 0.02% Sodium Azide, 2% BSA, PH 7.2
<b>Concentration</b>	0.5 mg/ml
<b>Storage</b>	Store at -20°C, Avoid freeze/thaw cycles
<b>Note</b>	For research use only.
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Affinity Purified
<b>Expiration Date</b>	12 months from date of receipt.