

Product Datasheet

PIK3R1 antibody orb258575

Description	Rabbit polyclonal antibody to PIK3R1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthetic peptide
Target	PIK3R1. This antibody has specificity with PIK3R1
Form/Appearance	PBS with 0.02% Sodium Azide, 2% BSA, PH 7.2
Concentration	0.5 mg/ml
Storage	Store at -20°C, Avoid freeze/thaw cycles
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Antigen Affinity Purified
NCBI	NP_852664
Expiration Date	12 months from date of receipt.