

Product Datasheet

Mouse IRF3 ELISA Kit orb437402

Description	Mouse Interferon Regulatory Factor 3 (IRF3) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10 ng/mL
Sensitivity	0.048 ng/mL
Tested Applications	ELISA
Target	IRF3
Storage	All reagents to be kept according to labels on the vials. TMB Substrate, Wash Buffer (30X concentrate), Stop Solution at 4°C; all others at -20°C.
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	tissue homogenates, cell lysates and other biological fluids.