

Product Datasheet

Substance P antibody orb303082

Description	Rabbit polyclonal antibody to Substance P
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IHC-Fr
Target	Substance P
Form/Appearance	PBS containing 15 mM sodium azide, pH 7.4
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, Avoid freeze/thaw cycles
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified
Expiration Date	12 months from date of receipt.