

Product Datasheet

Human ABCA3 ELISA Kit orb437640

Description	Human ATP Binding Cassette Transporter A3 (ABCA3) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Sensitivity	29pg/mL
Tested Applications	ELISA
Target	ABCA3
Storage	All reagents to be kept according to labels on the vials. TMB Substrate, Wash Buffer (30X concentrate), Stop Solution at 4°C; all others at -20°C.
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	tissue homogenates and other biological fluids.