

## Product Datasheet

### Human ADAM9 ELISA Kit orb437720

<b>Description</b>	Human A Disintegrin And Metalloprotease 9 (ADAM9) Sandwich ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156-10 ng/mL
<b>Sensitivity</b>	0.064 ng/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	ADAM9
<b>Storage</b>	All reagents to be kept according to labels on the vials. TMB Substrate, Wash Buffer (30X concentrate), Stop Solution at 4°C; all others at -20°C.
<b>Note</b>	For research use only.
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	tissue homogenates, cell lysates, cell culture supernates and other biological fluids.