

## Product Datasheet

### 14-3-3 antibody (HRP) (orb470128)

<b>Description</b>	Rabbit polyclonal antibody against 14-3-3 conjugated to HRP...
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Predicted Reactivity</b>	Mouse, Human
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	IHC-P, WB
<b>Immunogen</b>	KLH conjugated synthetic peptide derived from human 14-3-3 (201-245aa)
<b>Form/Appearance</b>	Lyophilized or Liquid
<b>Concentration</b>	- 100 µg (in 100 µl): 1 mg/ml
<b>Storage</b>	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
<b>Note</b>	For research use only.
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	27kDa
<b>Uniprot ID</b>	<a href="#">Q9CQV8</a> , <a href="#">P35213</a> , <a href="#">P68250</a> , <a href="#">P31946</a>
<b>Dilution Range</b>	WB: 1:50-200, IHC-P: 1:50-200
<b>Expiration Date</b>	12 months from date of receipt.