

Product Datasheet

Human ATG16L1 ELISA Kit orb438115

Description	Human Autophagy Related Protein 16 Like Protein 1 (ATG16L1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10 ng/mL
Sensitivity	0.063 ng/mL
Tested Applications	ELISA
Target	ATG16L1
Storage	All reagents to be kept according to labels on the vials. TMB Substrate, Wash Buffer (30X concentrate), Stop Solution at 4°C; all others at -20°C.
Note	For research use only.
Sample Types	tissue homogenates and other biological fluids.
Expiration Date	Please enquire.
Application Notes	tissue homogenates and other biological fluids.