

Product Datasheet

Rat CDNF ELISA Kit orb438554

Description	Rat Cerebral Dopamine Neurotrophic Factor (CDNF) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10 ng/mL
Sensitivity	0.067 ng/mL
Tested Applications	ELISA
Target	CDNF
Storage	All reagents to be kept according to labels on the vials. TMB Substrate, Wash Buffer (30X concentrate), Stop Solution at 4°C; all others at -20°C.
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	serum, plasma, tissue homogenates and other biological fluids.