

## Product Datasheet

### IVD antibody (HRP) orb477153

<b>Description</b>	Rabbit polyclonal antibody against IVD conjugated to HRP
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	ELISA, ICC, IHC-P, WB
<b>Immunogen</b>	KLH conjugated synthetic peptide derived from human IVD
<b>Target</b>	IVD (HRP)
<b>Form/Appearance</b>	Liquid: 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
<b>Concentration</b>	- 100 µg (in 100 µl): 1 mg/ml
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only.
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Affinity purified by Protein A
<b>MW</b>	43 kDa
<b>Uniprot ID</b>	<b>P26440</b> (Human); <b>Q9JHI5</b> (Mouse); <b>P12007</b> (Rat)
<b>Entrez</b>	<b>3712</b>
<b>Dilution Range</b>	WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200, ICC: 1:50-200
<b>Expiration Date</b>	12 months from date of receipt.