

## Product Datasheet

### Arfaptin 2 antibody orb10133

<b>Description</b>	Rabbit Polyclonal to Arfaptin 2.
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Predicted Reactivity</b>	Bovine, Canine, Sheep
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IF, IHC-P, WB
<b>Immunogen</b>	KLH conjugated synthetic peptide derived between 231-331 amino acids of human Arfaptin 2
<b>Target</b>	Arfaptin 2
<b>Form/Appearance</b>	Liquid: 0.01M TBS (pH7.4), 0.5% BSA, 0.015% ProClin™ 300 and 25% glycerol
<b>Concentration</b>	- 100 µg (in 200 µl): 0.5 mg/ml - 200 µg (in 400 µl): 0.5 mg/ml
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only.
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Affinity purified by Protein A
<b>MW</b>	38 kDa
<b>Uniprot ID</b>	<b>P53365</b> (Human); <b>Q8K221</b> (Mouse); <b>Q6AY65</b> (Rat)
<b>Entrez</b>	<b>23647</b>
<b>Dilution Range</b>	WB: 1:500-2000, ELISA: 1:500-1000, IHC-P: 1:400-800, IF: 1:100-500
<b>Expiration Date</b>	12 months from date of receipt.