

Product Datasheet

Goat IgG orb22143

Description	Goat IgG isotype control
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	DOT, IEP
Form/Appearance	Lyophilized, stable white powder. No preservative added.
Storage	This product is shipped lyophilized at ambient temperature. After arrival store at +4°C. Prolonged storage may be at - 20°C. Allow to stand at ambient temperature for 5-10 minutes to reach equilibrium. Reconstituted solution may be stored frozen. Working dilutions may be prepared by adding the required amount of phosphate buffered saline (PBS), should not be refrozen and preferably used within 24 hours. Repeated thawing and freezing should be avoided. If a slight precipitation occurs upon storage, this should be removed by centrifugation and will not affect the performance of the product.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.