

Product Datasheet

Trichomonas vaginalis antibody orb23579

Description	Trichomonas vaginalis is a mouse monoclonal antibody to infectious agent/toxin.
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, Purification
Immunogen	Disrupted Trichomonas vaginalis
Target	Trichomonas vaginalis. Trichomonas vaginalis
Form/Appearance	0.01M PBS, pH 7.2
Preservatives	0.1% Sodium azide
Storage	Short term up to 2 weeks store at -2 to -8°C. Long term, aliquot and store at to -20°C. Avoid multiple freeze/thaw cycles.
Note	For research use only.
Isotype	IgG3
Clonality	Monoclonal
Purity	Greater than 90% pure. Protein A chromatography
Source	Ascites
Biological Origin	Sample obtained from Mouse.
Clone ID	B989M
Expiration Date	12 months from date of receipt.
Application Notes	Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded. Capture: orb23579 Detection: orb23578