

Product Datasheet

Substance P antibody (FITC) orb53037

Description	Rabbit polyclonal antibody against Substance P conjugated to FITC.
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Tested Applications	ELISA
Immunogen	Recombinant Human Protachykinin-1 protein (1-129AA)
Target	Substance P (FITC)
Form/Appearance	Liquid: Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	>95%, Protein G purified
Uniprot ID	P20366
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.