

## Product Datasheet

### BRAF antibody orb394326

<b>Description</b>	Mouse monoclonal antibody to BRAF.
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC-P, WB
<b>Immunogen</b>	Recombinant fragment of human BRAF expressed in E. Coli
<b>Target</b>	BRAF. BRAF antibody detects endogenous levels of total BRAF
<b>Form/Appearance</b>	PBS (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Storage</b>	Store at -20°C
<b>Note</b>	For research use only.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Affinity-chromatography
<b>MW</b>	84.4 kDa
<b>Uniprot ID</b>	<b>P15056</b>
<b>Dilution Range</b>	ELISA: 1:10000, WB: 1:500-2000, IHC-P: 1:200-1000
<b>Expiration Date</b>	12 months from date of receipt.