

Product Datasheet

Rat IgG1,IgG2a,IgG2b, IgG2c heavy chain antibody (Biotin) orb414352

Description	Mouse monoclonal antibody against Rat IgG1,IgG2a,IgG2b, IgG2c heavy chain conjugated to Biotin
Species/Host	Mouse
Reactivity	Rat
Conjugation	Biotin
Tested Applications	ELISA
Target	Rat IgG1,IgG2a,IgG2b, IgG2c heavy chain (Biotin)
Form/Appearance	PBS: 0,1% NaN ₃ , 50% g lycerol
Storage	Store at -20°C.
Note	For research use only.
Isotype	IgG1, IgG2a kappa
Clonality	Monoclonal
Clone ID	MARG-COC-1
Expiration Date	12 months from date of receipt.