

## Product Datasheet

### Rat IgG1 isotype control (RPE) orb343824

<b>Description</b>	Rat IgG1 Isotype control conjugated to RPE
<b>Reactivity</b>	Rat
<b>Conjugation</b>	RPE
<b>Tested Applications</b>	ELISA, FACS, ICC, IF
<b>Immunogen</b>	Recombinant protein
<b>Target</b>	Rat IgG1
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1.0 mg/mL
<b>Preservatives</b>	0.01% (w/v) Sodium Azide
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Immunoaffinity chromatography
<b>Expiration Date</b>	12 months from date of receipt.