

Product Datasheet

Rat IgG2a isotype control (RPE) orb343825

Description	Rat IgG2a Isotype control conjugated to RPE
Reactivity	Rat
Conjugation	RPE
Tested Applications	ELISA, FACS, ICC, IF
Immunogen	Recombinant protein
Target	Rat IgG2a
Form/Appearance	Liquid
Concentration	1.0 mg/mL
Preservatives	0.01% (w/v) Sodium Azide
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG2a
Clonality	Monoclonal
Purity	Immunoaffinity chromatography
Expiration Date	12 months from date of receipt.