

Product Datasheet

ROBO-1 antibody orb345475

Description	Rabbit polyclonal antibody to ROBO-1
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Immunogen	Synthetic peptide (Internal)
Target	ROBO-1
Form/Appearance	Liquid: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Concentration	1.0 mg/mL
Preservatives	0.01% (w/v) Sodium Azide
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified
Uniprot ID	Q2M1J3
Entrez	6091
Expiration Date	12 months from date of receipt.
Application Notes	This affinity purified antibody has been tested for use in ELISA , IP, immunohistochemistry, immunofluorescence, and western blot. Specific conditions for reactivity should be optimized by the end user. Expect a band at ~181 kDa in size corresponding to ROBO-1 by western blotting in the appropriate cell lysate or extract.