

Product Datasheet

Animal IgG3 isotype control orb343833

Description	Animal IgG3 isotype control
Reactivity	Animal
Conjugation	Unconjugated
Tested Applications	ELISA, FACS, WB
Immunogen	Protein
Target	Animal IgG3
Form/Appearance	Liquid: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Concentration	1.0 mg/ml
Preservatives	0.01% (w/v) Sodium Azide
Storage	Store vial at -20° C prior to opening. Aliquot contents and freeze at -20° C or below for extended storage. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Note	For research use only.
Isotype	IgG3
Purity	Llama IgG3 isotype control has been prepared from llama serum by multiple chromatography steps using a combination of protein A and protein G chromatography. Coomassie stained SDS-PAGE of non-reduced llama IgG3 shows a band of ~100 kDa whereas the reduced
Expiration Date	12 months from date of receipt.