

Product Datasheet

Animal IgG1 isotype control (BAC) orb343841

Description	Animal IgG1 isotype control (BAC)
Reactivity	Animal
Conjugation	Unconjugated
Tested Applications	ELISA, FACS, WB
Immunogen	Protein
Target	Animal IgG1 (BAC)
Form/Appearance	Lyophilized: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free
Concentration	1.0 mg/mL
Preservatives	0.01% (w/v) Sodium Azide
Storage	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Note	For research use only.
Reconstitution	Reconstitution in 100 µl in Restore with deionized water (or equivalent)
Isotype	IgG1
Purity	Llama IgG1 isotype control has been prepared from llama serum by multiple chromatography steps using a combination of protein A and protein G chromatography. Coomassie stained SDS-PAGE of non-reduced llama IgG1 shows a band of ~150 kDa whereas the reduced
Expiration Date	12 months from date of receipt.