

Product Datasheet

Integrin beta 1 antibody orb348876

Description	Mouse monoclonal antibody to Integrin beta 1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS, IHC-Fr
Immunogen	Recombinant protein
Target	Integrin beta 1. Recognises human CD29.
Form/Appearance	PBS with preservative (0.02% Proclin 300)
Concentration	See vial label
Storage	Store at 4°C for up to 3 months. For longer term storage aliquot into small volumes and store at -20°C
Note	For research use only.
Isotype	Light chain kappa
Clonality	Monoclonal
Purity	Protein A affinity purified
Clone ID	K20
Expiration Date	12 months from date of receipt.
Application Notes	This antibody binds to CD29.