

## Product Datasheet

### XAGE1A antibody (Biotin) orb357604

<b>Description</b>	Rabbit polyclonal antibody against XAGE1A conjugated to Biotin
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	ELISA, IHC-P
<b>Immunogen</b>	Recombinant human X antigen family member 1 protein
<b>Target</b>	XAGE1A (Biotin)
<b>Form/Appearance</b>	Liquid: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only.
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	>95%, Protein G purified
<b>Uniprot ID</b>	<b>Q9HD64</b>
<b>Expiration Date</b>	12 months from date of receipt.